



PRODUCT
CATALOG



portacheck.com

1 Whittendale Drive Suite E, Moorestown, NJ 08057, info@portacheck.com, (856) 231-8894



The Health of the Herd is in Your Hands.



Since 2004, PortaCheck dairy animal testing products have made it easy and affordable to improve the health and production of herds around the world.

PortaCheck's reliable on-farm tests improve animal health and milk quality by detecting subclinical ketosis, mastitis, and other maladies before they become problems. PortaCheck is the recipient of three European innovative product development awards.

PortaCheck products are sold through our network of fine distributors. **To find a distributor near you, visit portacheck.com.**



BHBCheck Plus

Blood Ketone & Glucose Monitoring System
Bovine Blood

An accurate, affordable test that uses one drop of blood to provide quantitative results in 5 seconds (BHB) and 7 seconds (glucose).

This system is used to detect subclinical ketosis in bovine blood. Ketone strips sold in three varieties (10, 25, and 50 strips per vial).



PortaBHB

Milk Ketone Test
Bovine Milk

A semi-quantitative dipstick for the detection of subclinical ketosis. Simply dip the test strip into a vial of fresh milk or strip directly onto the test. Just one minute later, compare the color change to the color chart. Strips sold in two varieties (25 or 100 strips per vial). No refrigeration required.



PortaSCC for Cows

45-Min. Somatic Cell Count Test

Bovine Milk

This easy-to-use test provides semi-quantitative results for milk somatic cell count in just 45 minutes. The reaction causes the test strip to change color; users simply compare the result to a color chart. Compatible with our NEW PortaReader for quantitative results! (See back cover.)



PortaSCC Quick Test for Cows

5-Min. Somatic Cell Count Test

Bovine Milk

Short on time? The easy-to-use PortaSCC Quick Test provides semi-quantitative results for milk somatic cell count in just five minutes! The reaction with milk causes the test strip to change color; users simply compare the result to a color chart.



UdderCheck

Udder Health Test

Bovine Milk

UdderCheck measures Lactate Dehydrogenase (LDH), an enzyme present in milk when cells are damaged from udder infection. LDH is correlated to somatic cell count (SCC) and is not as easily affected by conditions such as stress, nutrition, parity, etc. In just two minutes, this milk dipstick test provides semi-quantitative results.



PortaSCC for Goats

Somatic Cell Count Test

Goat Milk

A somatic cell count test for early on-farm mastitis detection in goats. This 45-minute test helps dairy goat producers improve milk quality and production. Sold in two varieties (12 or 36 tests per kit). Results are read with a color chart.



PortaCheck products are sold through our network of fine distributors.

Find a distributor at portacheck.com

PortaReader

Clear, Quantitative Results for Mastitis Screening

Milk quality starts with healthy animals. The PortaReader interprets the results of PortaSCC somatic cell tests for subclinical mastitis, one of the leading causes of lost dairy profits. Increase detection accuracy with the PortaReader's clear quantitative results and easy-to-interpret display. This small, rugged unit can be carried in your pocket and used reliably in the field without breaking your budget.



Increase your milk production and profit: with the PortaReader, the health of your herd is in your hands!

**PortaCheck products are
sold through our network of
fine distributors.**

**To find a
distributor
near you, visit
portacheck.com.**

PortaCheck

1 Whittendale Drive
Suite E
Moorestown, NJ 08057 USA

info@portacheck.com
Tel 856.231.8894